

**Medical electrical equipment –
Part 2-44: Particular requirements for the safety
of X-ray equipment for computed tomography**

CORRIGENDUM 1

Page 22

29.1.105 Safety measures against excessive X-RADIATION

In the list of lettered items, instead of:

- d) When LOADING has been terminated by circumstances not noted under a), b) and c) above, resetting of the CT CONDITIONS OF OPERATION shall be required prior to the initiation of another scan.

read:

- d) When LOADING has been terminated by circumstances noted under a), b) and c) above, resetting of the CT CONDITIONS OF OPERATION shall be required prior to the initiation of another scan.